

User guide

VIKING®



Euro high back rehab chair

Product code 626S

Intended use

This chair is ideal for the dining room, lounge, bedroom etc where a standard chair is not suitable or the user has difficulty standing or sitting. Contoured backrest provides lumbar and neck support. Padded arms for extra comfort. Seat height, depth and angle are adjustable. This chair is suitable for domestic usage.



WARNING! This user manual must be read before using the chair. Severe injury or death may result if user instructions, maintenance instructions and product warnings are not followed.

Safe use

- **Do not exceed the safe working load of 200kg.**
- Ensure that the chair is stable on the floor and cannot slide.
- Inspect at regular intervals for wear and damage. Remove from use immediately if faults are found.
- Under no circumstances should the product be moved whilst occupied.
- Supervision or assistance when using the chair is recommended for users with limited mobility.
- The product is not suitable for use in bathrooms or wet areas.
- Only use this product indoors.
- Ensure there is adequate space around the chair to enable safe movement on, off and around the chair.
- Use the chair only for the purpose it is intended for.



WARNING! The decision to use this chair should always be guided by a qualified Health Care Professional, who can also advise and assist with the correct height adjustment and usage instruction.

Instructions for use

Your chair needs to be set to the required height under the guidance of a health care professional. This should enable you to sit and stand with ease.

Adjusting the seat depth.

To adjust the seat depth reach under both the left and right hand sides of the seat and unscrew each of the adjustable knobs. Push the back of the chair to the desired seat depth, then ensure you retighten both adjustable knobs under the seat.

Adjusting the seat height.


To adjust the seat height push the button in and slide the external leg piece up or down to the appropriate hole. Ensure the button protrudes by 2.5mm to 3mm to select the required height.



WARNING! Ensure all four legs are set at the same height. Severe injury or death may result if this is not done correctly

How to use the chair

All four legs must be level on the floor and cannot slide.



To sit on the chair, move backwards until you can feel the seat behind your legs, place each hand on the armrest. Support your body and prevent the chair from slipping away. Then lower yourself down.

When getting off the chair, move your bottom to the front edge of the seat ensuring your feet are firmly placed on the floor. Use the armrests to push yourself into a standing position leading with your head and shoulders.

Product specifications

Safe working load	200kg
Adjustable seat height	400mm - 550mm
Width inside arms	490mm
Seat width	500mm
Adjustable seat depth	480mm - 520mm
Height from seat to top of back	750mm

Options



Round base
glide feet
Code: 304-28FITTED



Serving tray
Code: 60T



Serving tray
storage bracket
Code: 60TS

Care and maintenance

- Clean and disinfect using regular household cleaning products.
- DO NOT use abrasive powders or scouring pads to avoid causing corrosion.
- If using Hyposal, care must be taken to remove all product residue with a clean soft cloth.
- Ensure chair is totally dry after cleaning.
- Should the equipment appear damaged, corroded or rusted remove from use immediately.

Disposal

Products that can no longer be used are to be disposed of separately from household rubbish. This must be done according to the local and national regulations for environmental protection and raw material recycling. Please look at your local council's website for further information on the correct disposal.

Expected lifetime

The expected lifetime of the chair is 3 years from the date of purchase.

Servicing

For any servicing, maintenance and troubleshooting queries, please contact the Cubro Service Department:

0800 656 527 +64 7 578 7228
sales@cubro.co.nz

149 Taurikura Dr, Tauriko, Tauranga, NZ

For additional copies of this user guide or for a different format, please contact Cubro.



WARNING! Do not make any changes or modifications to this product without consultation from Cubro.

Limited Warranty

1. This warranty applies only to the original purchaser of CUBRO LTD's product (who must be a purchaser who is resident or carrying on business in New Zealand and who has purchased the product directly from CUBRO LTD, or from a CUBRO LTD authorized distributor or reseller in New Zealand) and it is not transferable to any other person or entity.
2. This warranty replaces all conditions, warranties or guarantees that might be implied by law in favour of any person, which are excluded to the fullest extent permissible by law. Where the Consumer Guarantees Act would apply but the purchaser acquires or holds themselves out as acquiring any product for business purposes, the guarantees provided under the Consumer Guarantees Act are excluded in relation to that product.
3. CUBRO LTD warrants the product to be free from defects in materials and workmanship from date of purchase for a period of one year (12 months). This warranty does not cover any damage, defect, expense or loss of any kind caused by accident, misuse, abuse, neglect, negligence, alteration or modification (which includes the use of unauthorised parts or attachments), improper service, repair by other than authorized personnel or any defects not related to materials or workmanship. Wear of components in normal operation and failures resulting therefrom are excluded from this warranty.
4. CUBRO LTD has no obligation to the purchaser and is not obligated to honour all or any part of this warranty unless the following procedure is followed by the purchaser:
 - a. Before making a warranty claim, the purchaser should ensure that the product is defective by following standard 'trouble shooting' procedures and be able to attribute the fault to a defect in materials and/or workmanship of CUBRO LTD;
 - b. If CUBRO LTD determines that the product is not defective, or that the warranty claim is otherwise invalid, CUBRO LTD shall charge the purchaser a handling and servicing fee as well as any freights costs.
 - c. Upon receiving notice from the purchaser of an alleged defect in a product (which defect should be reported to CUBRO LTD immediately), CUBRO LTD will issue a serialized return management authorization (RMA). The purchaser must then return the entire unit or remove, at the purchaser's cost, the defective component part(s) identified, pack the unit or the component part(s) in a manner to avoid shipping damage and to ship the unit or the component part(s) to either CUBRO LTD, or a service centre as specified by CUBRO LTD, within 30 days of the date of the serialized return authorization date.
 - d. If CUBRO LTD require additional information relating to the use of the product by the purchaser or any other relevant information, the purchaser will supply such information as soon as practically possible and in such a manner as CUBRO LTD may reasonably require.
 - e. If CUBRO LTD access the purchaser's claim, it will either repair or replace the product, or repair or replace the defect in workmanship, as it may determine in its sole discretion.
5. Proper selection of a specific product for a specific application and operating environment, and its compatibility with other equipment is the purchaser's responsibility. CUBRO LTD does not warrant the performance of its products or their suitability for a particular purpose.
6. CUBRO LTD shall not be liable under any heading (whether in contract, negligence, or otherwise) for any consequential, indirect or incidental loss or damages of any kind and the maximum aggregate liability of CUBRO LTD (under any one or more headings) shall be the purchase price of the product.