



La-Z-Boy® Pinnacle platinum powerlift chair



Charcoal
Fabric



Twilight
Fabric



Flint
PU



La-Z-Boy® recliners are second to none when it comes to quality and style. The Pinnacle Platinum Power Lift Recliner makes life easier, and a lot more comfortable. Just press a button to get in or out, or fully recline all the way back to the zero-gravity position. The one-piece seat and back move together, eliminating gaps and providing full body support in all reclined positions. Features include power tilt and independent back and legrest motion. Plus, an easy-to-use wired controller and side storage pocket.

La-Z-Boy® chairs come with a special stain resistant iClean™ fabric, designed with innovative technology that creates a barrier of protection to repel spills before they turn into stains. Every inch of fabric is completely immersed in a stain resistant solution which guarantees every single fibre has been treated, leaving a concentrated barrier of protection but still keeps the fabric soft and luxurious.

Different fabrics and colours available.

Features

- Twin motor designed for lift and recline functions
- When fully reclined the user reaches zero-gravity position
- Back and legrest work together or independently
- Each chair function operates independently to achieve the best position for the resident
- Elevates the feet in any seating position not just in recline
- Easy-to-use wired controller and side storage pocket
- Full recline position stimulates blood circulation

Models	Colour	Codes
La-Z-Boy® Pinnacle fabric platinum powerlift chair	Charcoal	40050
La-Z-Boy® Pinnacle fabric platinum powerlift chair	Twilight	40051
La-Z-Boy® Pinnacle PU platinum powerlift chair	Flint	40060

Variety	La-Z-Boy® Pinnacle fabric platinum powerlift chair - Charcoal	La-Z-Boy® Pinnacle fabric platinum powerlift chair - Twilight	La-Z-Boy® Pinnacle PU platinum powerlift chair - Flint
Codes	40050	40051	40060
Warranty	Fabric: 3 years, Chair: limited lifetime	Fabric: 3 years, Chair: limited lifetime	Fabric: 3 years, Chair: limited lifetime
SWL	160kg	160kg	160kg
Overall Dimensions	W 830 D 990 H 1040 mm	W 830 D 990 H 1040 mm	W 830 D 990 H 1040 mm
Arm Height	N/A	N/A	N/A
Body Width	830 mm	830 mm	830 mm
Body Depth	990 mm	990 mm	990 mm
Body Height	1040 mm	1040 mm	1040 mm
Seat Width	510 mm	510 mm	510 mm
Seat Depth	550 mm	550 mm	550 mm
Seat Height	560 mm	560 mm	560 mm
Weight	66.7 kg	66.7 kg	66.7 kg



Made in America

Get in touch & connect   

Email sales@cubro.co.nz • Call us 0800 101 211 • Visit cubro.co.nz

Cubro®

© Copyright 2023 Cubro Ltd.