

Milano® overbed table

Code: 594

The Milano® overbed table features a one-piece, easy-to-clean top designed for superior infection control.

Suitable for use over both beds and chairs, it is easily height adjustable and its rounded edge helps to prevent items from falling off while the user is eating or reading.



FEATURES

- Designed to fit over most beds and chairs
- Seamless one-piece, moulded ABS plastic top for superior infection control
- Gas-assisted height adjustment that can be activated with one hand
- Soft rounded edges helps prevent items from slipping off
- 4x 50mm twin braking castors
- Easy to clean and made to last

Ready to order?

Our customer support team are here to help.

0800 656 527 | +64 7 578 7228 | sales@cubro.co.nz

Milano® overbed table

SPECIFICATIONS

Code(s)	594
Safe working load	20kg
Overall dimensions	Table top: W 430mm L 900mm
Height	725 - 1050mm
Inside base width	750mm
Castors	4x 50mm twin braking
Product weight	14.6kg
Material	ABS Plastic