

Roma® overbed pole

Code: 512

The free-standing overbed pole helps individuals to sit up and reposition themselves in bed. Made of chrome plated tubular steel with a moulded plastic handle, it is strong and durable, and can be flat-packed for easy storage. The easy to adjust strap gives a variation of approximately 600mm.

FEATURES

- Supports mobility and independence in bed
- Made of chrome plated tubular steel with chain strap and steel handle
- The easy to adjust chain gives a variation of approximately 600mm



Ready to order?

Our customer support team are available 5 days a week to help answer all your questions.

Phone: 0800 656 527

Phone: +64 7 578 7228

Email: sales@cubro.co.nz

Roma® overbed pole

SPECIFICATIONS

Code(s)	512
Safe working load	80kg
Overall dimensions	H 1820mm
Base dimensions	W 830 L 915mm
Handle height	Min 530 Max 680mm
Material	Chrome plated steel