

# Viking® shower chair with arms

Code: 334

The Viking® shower chair features an adjustable seat height and angle. It has been made from epoxy-coated mild steel and pre-treated for rust protection.



## FEATURES

- Made from epoxy-coated mild steel
- Drainage slots keep the seat free of water
- Ergonomically angled arms assist transfers
- Electrophoretic surface treatment to prevent rusting
- Non-slip rubber tips
- Height adjustable

**Ready to order?**

**Our customer support team are here to help.**

**0800 656 527 | +64 7 578 7228 | [sales@cupro.co.nz](mailto:sales@cupro.co.nz)**

# Viking® shower chair with arms

## SPECIFICATIONS

Code(s)	334
Overall dimensions	W 550mm D 590mm
Safe working load	125kg
Seat dimensions	W 440mm D 400mm
Seat height	400 - 550mm
Width between arms	450mm